

Reyoungel[®]



BIOHA

INTRODUCTION

- As a pioneer in medical aesthetics field, Bioha Laboratories is specialized in designing, manufacturing, and distributing implantable and injectable medical devices in the fields of plastic and aesthetic surgery, dermatology, and anti-aging.
- The SAX-HA (ZL 2012 1 0372786.0) patented technologies of Reyoungel® injectable cross-linked hyaluronic acid fillers have been developed by scientists and engineers specialized in biopolymers injectable medical devices with over 20-year full experience from Switzerland, France.
- Through our continuous research and development, we have developed the unique SAX technology patented in the United States.

When ageing, our face gradually develops fat and volume deficiency; our skin grows short of collagen and elasticity. Our products are designed to smooth outlines and wrinkles, reshape facial contours, and add volume to the lips. Our products are safe and effective when treating the entire face, creating a youthful, glamorous look for our customers.
- The quality management system of Bioha Laboratories is certified against ISO 13485 and Reyoungel® has obtained CE certification from the EU.
- Since the inception of Bioha Laboratories in 2003, millions of boxes of Reyoungel® products have been sold in over 70 countries worldwide, Reyoungel® also provides consumers with a safe, effective, and elegant aesthetic experience.

* Patent number ZL 2012 1 0372786.0
PCT No: PCT/CN2013/080217

BRAND HISTORY

Bioha Laboratories, established in 2003, is the mother company of Reyoungel, focusing on making dermal fillers.

Our facility is ISO 13485 registered and CE approved.

We are dedicated to the development of medical and beauty products. Our products range from Dermal Filler to PDO Thread, Mesotherapy, et cetera.

REYOUNGEL creates a series of dermal fillers used to treat signs

of ageing. They are composed of cross-linked hyaluronic acid extracted from

a non-animal- based HA gel, which is biodegradable and has a clear, colourless texture.



CONTENTS

1 **Fine Lines**

2 **Derm**

3 **Derm Deep**

4 **Derm Plus**

5 **Sub Skin**



R001

Fine Lines

1ml/2ml

Indications:

The correction of very thin wrinkles, eye-around wrinkles, "crow's feet" wrinkles around contour of mouth, etc.

PARAMETER

- Ingredient: Cross-linked hyaluronic acid (HA) gel
- HA Concentration: 20mg/ml
- Particle Size: 0.10-0.15mm
- Needle size: 30G/30G
- Duration: 6-9 months



NEEDLE
30G/30G





SUGGESTED CANNULA GAUGE
30G/30G



VOLUME
1ml/2ml

Injection depth: Upper part of dermis.





R002

Derm

1ml/2ml

Indications:

Correction of medium wrinkles, skin folds and the red border of the lips.

PARAMETER

- Ingredient: Cross-linked hyaluronic acid (HA) gel
- HA Concentration: 20mg/ml
- Particle Size: 0.15-0.28mm
- Needle size: 30G/30G
- Duration: 6-12 months



NEEDLE
30G/30G




SUGGESTED CANNULA GAUGE
30G/30G



VOLUME
1ml/2ml

Injection depth: Middle layer of the dermis.





R003

Derm Deep

1ml/2ml/10ml

Indications:

Medium and deep wrinkles filling, lips and cheek volume expansion, fill nasolabial folds.

PARAMETER

- Ingredient: Cross-linked hyaluronic acid (HA) gel
- HA Concentration: 20mg/ml
- Particle Size: 0.28-0.5mm
- Needle size: 26G/27G
- Duration: 6-12 months



NEEDLE

26G/27G



SUGGESTED CANNULA GAUGE

26G/27G



VOLUME

1ml/2ml/10ml

Injection depth: Deeper layer of the dermis.





R004 Derm Plus

10ml

Indications:

Derm plus is designed to fill deep wrinkles, recover and form volume of cheeks, cheekbones and chin, body shaping.

PARAMETER

- Ingredient: Cross-linked hyaluronic acid (HA) gel
- HA Concentration: 20mg/ml
- Particle size: 0.5-1.25mm
- Duration: 9-18 months
- Needle size: 23G/25G



NEEDLE
23G/25G



SUGGESTED CANNULA GAUGE
23G/25G



VOLUME
10ml

Injection depth: Inject with needle or cannula through an incision into the deeper layers of the dermis and/or subcutaneously into the skin depression that needs to be corrected.





R005 Sub Skin

2ml/10ml

Indications:

Sub Skin is designed to form volume, correct facial contours, expand volume of cheeks, chin, restore and change of nasal contour, is used to compensate for the loss of volume during skin aging, body shaping.

PARAMETER

- Ingredient: Cross-linked hyaluronic acid (HA) gel
- HA Concentration: 20mg/ml
- Particle Size: 1.25-2.0mm
- Needle size: 22G/23G
- Duration: 9-18 months



NEEDLE
22G/23G




SUGGESTED CANNULA GAUGE
22G/23G



VOLUME
10ml

Injection depth: Inject with cannula through the incision or with needle 22G/23G Subcutaneous or periosteal into the area of skin depression that needs to be corrected.



CERTIFICATES

European CE Mark CE2764

EC Design-Examination Certificate

EN ISO 13485: 2016



GLOBAL SALES NETWORK

Specialized In The Field Of Sodium Hyaluronate

Acid Gel Around 20 Years



ADVANTAGES (Technology)

CROSS-LINKING TECHNOLOGY

SAX-HA" technology is a new type of cross-linking technology that enables synthetic hydrogels with hydrophobicity on the polymer chain to obtain strong self-association functions through strong hydrophobic association. This self-association is dynamic, and injected gels tend to stay and clump together to maintain shape and reduce gel migration in the body.

There are two commonly used crosslinking agents:


BDDE (butylene glycol glycidyl ether) and DVS (divinyl sulfone).

They have certain residues and adverse reaction reports.

Reyoungel uses **SAX-HA** cross-linking technology, which has the following advantages:

1. More thorough removal of cross-linking agent (**BDDE**) residues, reducing adverse reactions.
2. Better retention, reducing displacement. Make the shaping effect more natural and longer-lasting(6-18 months).

Cross-Linking - HCXL



STUDY DATA


A pre-market clinical trial of Reyoungel followed up 357 patients from China and Russia and found that due to Reyoungel's lower in vitro enzymatic hydrolysis rate, its effect is more lasting than similar products in Europe.

Clinical Data




Reyoungel completed 357 cases of clinical trials before marketing. All of these cases spread in China and Russia.

In order to get clinical data continually. Reyoungel built a well-established post market clinic follow-up systems.



Reyoungel VS other EU brand
Enzymolysis in vitro



Reyoungel enzymolysis in vitro time is less than other similar brands, this data shows that Reyoungel has longer lasting.

Reyoungel®



UK HEADQUARTERS

BIOHA LABORATORIES GROUP LTD

9 PANTYGRAIGWEN ROAD PONTYPRIDD

MID GLAMORGAN ENGLAND

UNITED KINGDOM CF37 2RR

Web: www.reyoungel.com

E-mail: info@reyoungel.com

